

# HLT37015 Certificate III in Sterilisation Services *outline*

There are 14 units of competency, including 10 core and 4 elective units required to attain qualification.

Code	Title	Type
HLTWHS001	Participate in workplace health and safety	Core
HLTWHS005	Conduct manual tasks safely	Core
HLTINF006	Apply basic principles and practices of infection prevention and control	Core
HLTSTE001	Clean and disinfect reusable medical devices	Core
HLTSTE002	Inspect and pack reusable medical devices	Core
HLTSTE003	Sterilise loads	Core
HLTSTE004	Manage sterile stock	Core
HLTSTE005	Care for reusable medical devices	Core
CHCCOM005	Communicate and work in health or community services	Core
CHCDIV001	Work with diverse people	Core
BSBFLM312	Contribute to team effectiveness	Elective
BSBFLM303	Contribute to effective workplace relationships	Elective
BSBWOR301	Organise personal work priorities and development	Elective
PSPGEN010	Contribute to workgroup activities	Elective

The course is self-paced over a maximum of 52 weeks duration. Depending on your experience and workplace flexibility, it may take between 24 and 52 weeks to complete. The course is comprised of online learning and assessment, as well as logbooks and workplace skills assessment activities. Workplace logbook activities must be completed under the direct supervision and observation from a suitably qualified colleague within the workplace.

To successfully complete this course, you will need to be able to set aside a minimum of 3 hours per week to complete online learning activities. To complete all practical activities for assessment, you will need to coordinate rotation through the different workstations in your sterilising service or reprocessing department to ensure you gain the skills, knowledge, and experience on all assessable equipment. Your final workplace skills assessment will be undertaken in the workplace.